

DETACHABLE ENDO RETRIEVAL

EASY TO USE, ONE-HANDED OPERATION

- Minimizes Contamination to the Body
- Contains Specimens and Stones
- Easy to Use Closing Mechanism
- Prevents Serious Infection
- Easy To Read Indicator for Bag Orientation
- Insertable Into a 10mm or Longer Canula
- Push/Pull String Loaded, Rod Activated Bag
- Not Made With Natural Rubber Latex



AVAILABLE IN THREE SIZES

TM2001	85mm x 135mm	(3"x6")	250 ml
TM2000	125mm x 135mm	(5"x7")	500 ml
TM2002	190.5mm x 228.6mm	(7.5" x 9")	1200 ml