

SURE-FLO INSUFFLATION TUBING

COMPATIBLE WITH MOST INSUFFLATORS

- 1/2" Tubing
- 10' Soft High-Flow Tubing
- 0.1 Micron Filter (99.99% Efficient)
- Hydrophobic ULPA Filter
- DEHP Free
- No Natural Rubber Latex



TM6500

Heated Insufflation Tubing

INSUFFLATION TUBING

CROSS REFERENCE

Description	Manufacturer	Item #	QTY per UOM	Description	Part #	QTY per UOM
High Flow Insufflation Tubing, Rotating Luer Lock, 0.3 Micron Filter	Carefusion / BD V.Mueller	LA103-RLL	Bx20	Insufflation Tubing 10ft, .01 micron filter with dual swivel luer connectors	TMITF10	Bx10
Insufflation Tubing (10ft length) & 0.1 micron hydrophobic filter assembly with luer lock connector.	ConMed	60-6050-099	Bx10	Insufflation Tubing 10ft, .01 micron filter with dual swivel luer connectors	TMITF10	Bx10
Insufflation Tubing (10ft length) & 0.1 micron hydrophobic filter assembly with suction type connector.	ConMed	60-6050-097	Bx10	Insufflation Tubing 10ft, .01 micron filter with swivel luer & funnel connector	TM6000	Bx10
Insufflation Tubing (10ft length) & 0.1 micron hydrophobic filter assembly with "O" ring connector.	ConMed	60-6050-098	Bx10	Insufflation Tubing 10ft, .01 micron filter with swivel luer & quck connector	TM102	Bx10
10' Clear Insufflation Tubing, 0.1 Micron Filter	DeRoyal	28-0207	Cs20	Insufflation Tubing 10ft, .01 micron filter with dual swivel luer connectors	TMITF10	Bx10
10' High Flow Insufflation Tubing, .3 Micron Filter with Stryker Funnel Connector	DeRoyal	28-0211	Cs20	Insufflation Tubing 10ft, .01 micron filter with swivel luer & funnel connector	TM6000	Bx10
10' Clear Insufflation Tubing w/Quick Connect Fitting, .1 Micron Filter	DeRoyal	28-0209	Cs20	Insufflation Tubing 10ft, .01 micron filter with swivel luer & quick connector	TM102	Bx10
10' High Flow Insufflation Tubing, .3 Micron Filter w/ Storz Fitting	DeRoyal	28-0212	Cs20	Insufflation Tubing 10ft, .01 micron filter with swivel luer & 15/22mm connector	TM1522	Bx10
10' High Flow Insufflation Tubing, .3 Micron Filter w/ Linvatec Fitting	DeRoyal	28-0216	Cs20	Insufflation Tubing 10ft, .01 micron filter with swivel luer & 15/22mm connector	TM1522	Bx10
Insufflation Tubing, 10ft, Luer, 0.1 Micron Filter	MediChoice Owens & Minor	3486INSF1001	Cs60	Insufflation Tubing 10ft, .01 micron filter with dual swivel luer connectors	TMITF10	Bx10
Insufflation Tubing, 10ft, Quick Connect 0.1 Micron Filter	MediChoice Owens & Minor	3486INSF1002	Cs60	Insufflation Tubing 10ft, .01 micron filter with swivel luer & quick connector	TM102	Bx10
Insufflation Tubing, 10ft, Funnel Connect 0.1 Micron Filter, Stryker	MediChoice Owens & Minor	3486INSF1003	Cs60	Insufflation Tubing 10ft, .01 micron filter with swivel luer & funnel connector	TM6000	Bx10
Insufflation Tubing, 10ft, Luer, Round Connect, 0.1 Micron Filter, Storz	MediChoice Owens & Minor	3486INSF1004	Cs60	Insufflation Tubing 10ft, .01 micron filter with swivel luer & 15/22mm connector	TM1522	Bx10
Dual Rotating Luer Connectors	Medline	DYNJ011933	Cs10	Insufflation Tubing 10ft, .01 micron filter with dual swivel luer connectors	TMITF10	Bx10
CPC Connector (Quick Connect) & Rotating Luer	Medline	DYNJ012933	Cs10	Insufflation Tubing 10ft, .01 micron filter with swivel luer & quick connector	TM102	Bx10
Push-on Connector & Rotating Luer	Medline	DYNJ014933	Cs10	Insufflation Tubing 10ft, .01 micron filter with swivel luer & funnel connector	TM6000	Bx10
22mm Round Connector and Rotating Luer	Medline	DYNJ016933	Cs10	Insufflation Tubing 10ft, .01 micron filter with swivel luer & 15/22mm connector	TM1522	Bx10
CPC Connector (Quick Connect) & Rotating Luer	Smith & Nephew	7205861	Bx10	Insufflation Tubing 10ft, .01 micron filter with swivel luer & quick connector	TM102	Bx10



REQUEST A SAMPLE
FROM YOUR REP