

Humeca Mesher and V-Carriers



Humeca Mesher - features

- Lightweight anodized aluminum machine
- Integrated blades and axle
- 50 parallel circular blades
- Spring construction to assure pressure on skin during meshing
- 2 configurations available; V10 and V15
- V10 version compatible with Zimmer® Dermacarrier II
- V15 version compatible with Aesculap®/B. Braun® carrier
- Continuous rotational drive
- Cutting axle can easily be replaced
- Sterilization case from stainless steel with silicone parts available designed for the Humeca Mesher for cleaning and sterilization
- Dimensions: l x w x h: 220x212x183 mm (8.7x8.3x7.2")
- Weight: 4.4 kg (9.7 lb)

V-carriers - features

- Available V-carrier ratios:
V10: 1:1, 1:1.5, 1:2, 1:3 and 1:6
V15: 1:1, 1:1.5 and 1:3
- V-pattern in all carriers (except 1:6). Groove patterns of V-carriers connect to each other.
- Flexible transparent polypropylene material, medical grade
- Individually sterile packed in peel pouch.
- Dimensions (l x w x h):
V10: 280 x 79 x 1.0 mm (11.0x3.1x0.04")
V15: 280 x 79 x 1.5 mm (11.0x3.1x0.06")

Order information

- | | |
|-----------|--|
| 6.V10-1.0 | V-carrier, perforation 1:1, for Humeca and Zimmer® meshers, box 10 pcs. |
| 6.V10-1.5 | V-carrier, expansion ratio 1:1.5, for Humeca and Zimmer® meshers, box 10 pcs. |
| 6.V10-2.0 | V-carrier, expansion ratio 1:2, for Humeca and Zimmer® meshers, box 10 pcs. |
| 6.V10-3.0 | V-carrier, expansion ratio 1:3, for Humeca and Zimmer® meshers, box 10 pcs. |
| 6.V10-6.0 | V-carrier, expansion ratio 1:6, for Humeca and Zimmer® meshers, box 10 pcs. |
| 6.V15-1.0 | V-carrier, perforation 1:1, for Humeca and Aesculap® / B.Braun® meshers, box 10 pcs. |
| 6.V15-1.5 | V-carrier, expansion ratio 1:1.5, for Humeca and Aesculap® / B.Braun® meshers, box 10 pcs. |
| 6.V15-3.0 | V-carrier, expansion ratio 1:3, for Humeca and Aesculap® / B.Braun® meshers, box 10 pcs. |
| 6.HM01 | Humeca® mesher |
| 6.213 | Mesher cutting axle with blade block |
| 6.HMAC03 | Autoclave case for Humeca® mesher and cutting axle 250x250x185 mm (9.8x9.8x7.3") |

