



Endotracheal Tube Electrodes

The Surface Electrode is designed to be secured to an endotracheal tube and to be positioned in contact with the true vocal cords to assist in monitoring the integrity of the Recurrent Laryngeal Nerve branch of the Xth Cranial Nerve during surgical procedures of the neck.

Supports virtually all nerve monitors

No additional cables needed

Universal ET tube compatibility

Superior electrical conductivity



Manufactured in the USA.

Endotracheal Tube Electrodes

Carotid Endarterectomy - Craniotomy - Thyroidectomy - Anterior Cervical Fusion - Skull Base Procedures

ORDERING INFORMATION

Thyroidectomy 2-Channel Kit



Catalog No.	Description	Quantity
29-80500DCK	One (1) each 29-80500DC-Dual Channel Laryngeal Electrode	5 kits/box
	One (1) each NE-S-250013/0.4- White Needle Electrode	
	One (1) each 29-80250- Green return needle electrode	

Thyroidectomy Laryngeal Kit



Catalog No.	Description	Quantity
29-80500K	One (1) each 29-80500 S- Standard Single Channel Electrode	5 kits/box
	One (1) each PNT 0.8/90 - Flush Tip Monopolar Stimulator Probe	
	One (1) each SEAg-S- 10000/20 Green Hydrogel Surface Electrode	

Stim Kit



Catalog No.	Description	Quantity
29-80500SK	One (1) each 29-80500 S- Standard Single Channel Electrode	5 kits/box
	One (1) each NE-S-250013/0.4 - White Needle Electrode	
	One (1) each - Red Disposable Stim Wire	
	One (1) each SEAg-S-10000/20 - Green Hydrogel Surface Electrode	

Dual Channel Electrodes



Catalog No.	Description	Quantity
29-80500DC	Standard Dual Channel Electrode, sterile, single use Fits 6.0 to 7.5 ET Tube	5

Single Channel Electrodes



Catalog No.	Description	Quantity
29-80500S	Standard Single Channel Surface Electrode, sterile, single use, Fits 7.5 to 10.0 ET Tube	5
29-80500W	Wrap Around Single Channel Surface Electrode, sterile, single use, Fits 7.5mm to 10mm ET Tube	5

Disposable Subdermal Needle Electrodes



Catalog No.	Description	Size	Quantity
NE-S-2500-13/0.4	White Single Subdermal Needle Electrode	13 x 0.40mm (0.5" x 27G)	5/Box
29-80250	Green Single Ground Needle Electrode	27 Gauge	5/Box
NE-T-1500/13/0.4R	Red/Red Single Twisted Subdermal Needle Electrode	13 x 0.40mm (0.5" x 27G)	5/Box
NE-T-1500/13/0.4B	Blue/Blue Single Twisted Subdermal Needle Electrode	13 x 0.40mm (0.5" x 27G)	5/Box

Probes



Catalog No.	Description	Quantity
PNT 0.8/90	Flush Tip Monopolar Stimulator Probe	10
20-3001	Disposable Monopolar Direct Nerve Stimulator Probe	10