

Our Sigmoidoscope Kits save time and eliminate the risk of cross contamination. Used in sigmoidoscopy procedures for visual through the sigmoid colon for locating ulcers, abnormal cells, polyps, and cancer.

Our sigmoidoscope and anoscope kits are single-use, 100% latex- free, require no filter, and eliminate the expense associated and complexity of cleaning.


- Patient safe
- Cost effective
- Fewer components
- LED light lasts 20+ hours
- Minimizes patient discomfort
- Designed to minimize tissue damage
- Eliminates the risk for cross contamination

Sigmoidoscope Kits

Products	Bioseal lighted Sigmoidoscope	W-A Sigmoidoscope	Bulb and Bladder	Sigmoidoscope suction with 6' connecting tube 18 Fr.	Procto Swabs
ASDKS01	*		*		1
ASDSK02	*		*	*	4



ASDKS01



ASDSK02

Surgical Mini Lighted Anoscope

ASDSA03

Size	Lighted Bevelled 10cm
UOM	1/pk 10pks/Cs
Sterile/Non-Sterile	Sterile
Latex/Latex-free	Latex Free
Instrument Type	Surgical Mini Lighted Anoscope
Disposable/ Single use	Disposable/ Single use
Department	Gastroenterology



Surgical Maxi Lighted Anoscope

ASDSA04

Size	(Lighted)
UOM	10 Ea/Cs
Sterile or Non-Sterile	Sterile
Latex or Latex-Free	Latex Free
Kit Components	1 Lighted Scope Recto Maxi 25cm
Instrument Type	Anoscope
Disposable or Reusable	Disposable/ Single use
Department	Gastroenterology



Lighted Surgical Anoscope

ASDLS05

Size	10cm Mini
UOM	1/pk 10pks/Cs
Sterile/Non-Sterile	Sterile
Latex/Latex-free	Latex Free
Instrument Type	Mini Lighted Surgical Anoscope
Disposable/ Single use	Disposable/single use
Department	Gastroenterology

